

The equine clinic of Conques saves more than 25,000 € per year with EtO gas sterilisation.



“With the Anprolene, I can reuse my material several times.”

Dr Serge Lenormand, Equine clinic of Conques,  
*Saint-Aubin de Branne, France*

Today, the Conque stud farm is and insemination centre with a world class reputation thanks to their cutting edge technology and high quality supporting infrastructure.

The Conque stud farm specialises in: breeding, insemination, embryo transfer and foaling.



Items :	Quantity (per month)	€ / unit	Good for how many autoclave cycles ?	Saving per month
EZ Way Filter	30	14.30	/	429 €
« Y » infusion drip	30	14.30	/	429 €
GIA stapler	4	290	/	1160 €
Laparoscopic electrodes	10	12	/	120 €
Others				200 € est.
2 cycles per month at 23€ per cycle give a final saving of 2,292 € per month.				

The Andersen sterilisation system releases gas directly to the heart of the load using its unique micro-dose technology. This reduces the amount of gas required and speeds up the process.

For more informations, contact :  
**HW Andersen Products**  
IICK  
Etienne Sabbelaan 51  
8500 Kortrijk, Belgium  
belgium@anderseneurope.com

SAVE WITH ANPROLENE®

AN74

Anprolene®  
gas sterilisation

AN74i with abator



### Design

The AN74i's compact table-top size (net weight 20.8kg) makes it suitable for a variety of small spaces. The cabinet is made of hard-wearing stainless steel and front loads for operator convenience.

The AN74ix has a longer cabinet to allow for the sterilisation of longer/larger items. The AN74ix cabinet has an interior length of 74.3cm, compared to 49.5cm for the AN74i.

### Installation

The AN74i can be easily installed either on a storage cabinet base, or simply placed upon a work surface. High pressure, low volume pumps simplify installation, requiring only a 7/8" reinforced PVC exhaust hose to the outside. No water or vacuum line required. Wall mount installation kit and hoses included.

### Safe for the Environment

EPA registered (Environmental Protection Agency), a combination of purge and ventilation systems ensures that EtO emission levels are kept well below prescribed levels. For additional care we offer the AN5100 Gas Abator which removes more than 98% of EtO emitted from the AN74.

### Certification

The AN74 is manufactured by an ISO certified company and is CE marked. CSA International tested for electrical safety.

### Warranty

This steriliser comes with a one-year warranty which, for an additional sum, may be extended to five years.

The customer may choose to pass the exhaust gases through an abator. This may be required by a strict local regulation, or the customer may be venting the exhaust into their central waste gas system rather than directly to the outside.

The unit uses a disposable cartridge, which will abate 200 17.6 gram ampoules. Like the gas ampoules, the abator cartridge is free of EtO after use and may be disposed of as household waste.

**Anprolene sterilisation cycles are designed to be the gentlest available, allowing repeated sterilisation whilst preserving medical devices in pristine condition.**

AN74ix



ANPRO®

ANDERSEN  
PRODUCTS