

100% Ethylene Oxide Sterilisation

Most effective sterilant  
Most efficient steriliser

AN85 Exposure Indicators  
AN86 Exposure Indicator Tape

**Anprolene®**  
gas sterilisation

**EOGas®**

**Indicators Stickers**  
Quality Assurance Products  
designed for use with  
Anprolene and EOGas sterilisers.

### AN85 Exposure Indicator Stickers

These labels provide immediate visual assurance that items have been exposed to Ethylene Oxide gas by undergoing a colour change from yellow to blue. They are suitable for use with Seal and Peel®, Tyvek®, CSR wrap and conventional sterilisation paper. Featuring strong adhesive backing they are attached to each item, before wrapping and sterilising.

Conforms to ISO 11140

**Size:** 2.5cm (W) x 4.7cm (L)

**Quantity:** 200/box

**Shelf life:** 2 years from date of manufacture



### AN86 Exposure Indicator Tape

This package sealing tape provides a reliable visual guarantee that items have been exposed to Ethylene Oxide gas by changing from yellow to blue. It is suitable for use with Seal and Peel®, Tyvek®, CSR wrap and conventional sterilisation paper. Featuring strong adhesive backing they are attached to each item, before wrapping and sterilising. There is a blank area to record information such as sterilisation date or lot number.

Conforms to ISO 11140

**Size:** 1.9cm (W) x 11.1cm (L)

**Quantity:** 200/box

**Shelf life:** 2 years from date of manufacture



*These products do not provide a sterility guarantee. They are not intended to replace Dosimeters (AN87, AN1087) or biological controls which measure the kill of bacterial spores (AN2203 and AN2202).*

Anprolene® and EOGas are registered trademarks of H.W. Andersen Products, Inc.

**ANPRO** **ANDERSEN**  
**PRODUCTS**  
the future of gas sterilisation

H.W. Andersen Products Ltd  
Unit 808, Oakwood Business Park  
Fowler Road  
Clacton-On-Sea, CO15 4AA  
United Kingdom  
Tel.: +44 (0)1255 428 328  
email: [uk@anderseneurope.com](mailto:uk@anderseneurope.com)  
[www.anderseneurope.com](http://www.anderseneurope.com)

H.W. Andersen Products, Inc.  
Health Science Park  
3202 Caroline Drive  
Haw River, NC 27258-9564  
USA  
Tel.: 800-523-1276  
email: [customerservice@anpro.com](mailto:customerservice@anpro.com)  
Ref: AN85/AN86 Flyer 201602a